




Summit Hill Laboratories  
3 Sheila Drive  
Tinton Falls, NJ 07724  
[sales@summithilllaboratories.com](mailto:sales@summithilllaboratories.com)

# VETROSON® MILLENNIUM MOTOR PACK

**A POWERFUL UNIT WHICH MAY BE USED AFTER SCALING TEETH TO CLEAN AND POLISH. MAY BE USED FOR SECTIONING AND RESHAPING TEETH. SELECTION OF BURS ENABLES VETERINARIAN TO PERFORM ANY ENDODONTAL WORK REQUIRED.**




**THE VETROSON® MOTOR PACK CONSISTS OF:**

- A COMPACT CONTROL UNIT WHICH DELIVERS MORE POWER THAN NECESSARY FOR POLISHING TEETH
- A LIGHT-WEIGHT 5" 20,000 RPM VARIABLE SPEED MICRO MOTOR HANDPIECE
- A SNAP-ON PROPHY ANGLE
- VETROSON® SNAP-ON POLISHING CUPS
- VETROSON® CLEANER/LUBRICANT SOLUTION

## **THE VETROSON® TRAUMA KIT IS AVAILABLE AS AN ACCESSORY**

The Trauma Kit consists of a latch type contra angle for holding the #2, #4, #6 & #8 round burs for cleaning out existing cavities, a #33 ½ & #35 inverted cone for shaping bur-cleaned cavities, a diamond cylinder for reshaping broken teeth, an 18mm x .45 single-sided diamond disc for cutting or sectioning, and a bur block and bur tray for cleaning burs.



**THE VETROSON® MOTOR PACK HAS BEEN DESIGNED TO BE USED IN CONJUNCTION WITH THE VETROSON® MILLENNIUM 25/30 ULTRASONIC SCALER.**

**FIVE YEAR WARRANTY**